



Integrating the Healthcare Enterprise (IHE) Integration Statement

Vendor	Product Name	Version	Date
Greenway Medical Technologies	PrimeSuite	PrimeSuite 2008	1 Feb 2010

This product implements all transactions required in the **IHE Technical Framework** to support the **IHE Integration Profiles, Actors** and **Options** listed below:

Integration Profiles Implemented		Actors Implemented	Options Implemented
APS	Antepartum Summary	Content Creator	None
		Content Consumer	View Option Document Import Option
CRD	Clinical Research Data Capture	Form Filler	None
EDR	Emergency Department Referral	Content Consumer	View Option Document Import Option
IC	Immunization Content	Content Creator	None
		Content Consumer	View Option Document Import Option
RFD	Retrieve Form for Data Capture	Form Filler	SOAP Transport Option Archive Form Option Data Clarifications Option
SVS	Sharing Value Sets	Value Set Consumer	HTTP Binding Option SOAP Binding Option
XPHR	Exchange of Personal Health Record Content	Content Creator	None
		Content Consumer	View Option Document Import Option
CT	Consistent Time	Time Client	None
DSC	Drug Safety Content	Form Filler	None
LDS	Labor and Delivery Summary	Content Consumer	View Option Document Import Option
PDQ	Patient Demographics Query	Patient Demographics Consumer	None
RPE	Retrieve Protocol for Execution	Protocol Executor	None
XDS-MS	Cross-Enterprise Sharing of Medical Summaries	Content Creator	Referral Option
		Content Consumer	View Option Document Import

Internet address for vendor's IHE information:

www.greenwaymedical.com/solutions/prime-suite/

Links to General Information on IHE

North America:	http://www.ihe.net
Europe:	http://www.ihe-europe.org
Japan:	http://www.jira-net.or.jp/ihe-j



Healthcare Information Technology Standard Panel (HITSP) Integration Statement

Vendor	Product Name	Version	Date
Greenway Medical Technologies	PrimeSuite	PrimeSuite 2008	1 Feb 2010

This product implements all transactions required by HITSP to support the **Interoperability Specifications, Transaction Packages, Transactions and Components** listed below:

Specification Implemented	Actors Implemented	Options Implemented
T-16 Consistent Time Transaction	Time Client	None
T-23 Patient Demographics Query Transaction	Patient Demographics Consumer	None
T-66 Retrieve Value Set Transaction	Value Set Consumer	HTTP Binding Option SOAP Binding Option
TP-50 Retrieve Form for Data Capture Transaction Package	Form Filler	SOAP Transport Option Archive Form Option Data Clarifications Option
C-32 Summary Documents Using HL7 Continuity of Care Document (CCD) Component	Content Creator	None
	Content Consumer	View Option Document Import Option
C-38 Patient Level Quality Data Document Using IHE Medical Summary (XDS-MS) Component	Content Creator	Referral Option
	Content Consumer	View Option Document Import
C-48 Encounter Document Using IHE Medical Summary (XDS-MS) Component	Content Creator	None
	Content Consumer	View Option Document Import
C-78 Immunization Document Component	Content Creator	None
	Content Consumer	View Option Document Import Option

Internet address for vendor's IHE information:

www.greenwaymedical.com/solutions/prime-suite/

Links to General Information on IHE

North America:	http://www.ihe.net
Europe:	http://www.ihe-europe.org
Japan:	http://www.jira-net.or.jp/ihe-j